



## Alt-R™ CRISPR-Cas9 System

Delivery of ribonucleoprotein complexes into HEK-293 cells using the Lonza® Nucleofector™ System

> SEE WHAT MORE WE CAN DO FOR YOU AT [WWW.IDTDNA.COM](http://WWW.IDTDNA.COM).

custom oligos • qPCR • next generation sequencing • RNAi • genes & gene fragments • CRISPR genome editing

# REVISION HISTORY

---

Version	Date released	Description of changes
4	June 2021	Updated to include additional information on new products.
3.1	August 2019	Specified unit of measure for centrifuge from 600 rpm to 150 x g.
3	July 2018	Added instructions for using Alt-R CRISPR-Cas9 sgRNA. Updated names and catalog numbers for Alt-R enzymes (V3). Updated to current IDT styles and formatting.
2.3	May 2018	Added note about use of improved Alt-R enzymes (V3): direct substitution in protocol of V3 enzymes for original enzymes (3NLS).
2.2	October 2017	Added information about new IDT crRNA design tools.
2.1	August 2017	Added information about new IDT products (Alt-R Cas9 variants).
2	January 2017	Updated product names to specify CRISPR-Cas9 system to differentiate these from CRISPR-Cpf1 system reagents. Replaced custom Ultramer oligo with the Alt-R Cas9 Electroporation Enhancer. Added ordering information and references for the fluorescently labeled tracrRNA, Alt-R CRISPR-Cas9 tracrRNA – ATTO 550. Corrected typographical errors.
1	November 2016	Original protocol.

# TABLE OF CONTENTS

Revision history	2
Introduction	4
Important considerations	4
Consumables and equipment	6
Equipment	6
Consumables from other suppliers	6
Consumables from IDT—Reagents	7
Protocol	8
Culture cells	8
Prepare RNA	8
Form the RNP complex	9
Prepare the Nucleofector™ system	9
Perform electroporation of cells with the Nucleofector™ system	9
References	11

# INTRODUCTION


This protocol describes the delivery of a CRISPR-Cas9 ribonucleoprotein (RNP) complex, containing Alt-R CRISPR-Cas9 guide RNA (crRNA:tracrRNA duplex or sgRNA) and a Cas9 enzyme (nuclease or nickase), into HEK-293 cells using electroporation with the Nucleofector™ system (Lonza) and is based on 2 protocols: the [Lonza 96-well Shuttle Protocol for HEK-293](#) [2] and the [IDT Alt-R CRISPR-Cas9 System User Guide](#) [3].

Go to [www.idtdna.com/CRISPR-Cas9](http://www.idtdna.com/CRISPR-Cas9) (Resources section, Application notes), for tips on using the following:

- Fluorescently-labeled tracrRNA (Alt-R Cas9 tracrRNA – 5'ATTO™ 550, 488, and 647, as well as GFP-Cas9) to investigate electroporation efficiency, or to select for transfected cells via cell sorting.
- Nickases (Alt-R *S.p.* Cas9 D10A Nickase or Alt-R *S.p.* Cas9 H840A Nickase) to help reduce off-target effects and promote homology-directed repair.

## Important considerations

1. Use low-passage, healthy cells. A critical factor affecting the success of electroporation is the health of the cells. It is important to:
  - Use the lowest passage number cells available
  - Subculture cells for at least 2–3 days before the electroporation procedure
  - Replace the media the day before electroporation
  - Determine the optimal confluency for your cell type

 **Note:** Optimal confluency for HEK-293 cells is 70–85% at the time of nucleofection. Higher cell densities are likely to reduce viability and electroporation efficiency [1].
2. Wash cells after trypsinization. Trypsin and FBS may contain RNase activity that can quickly degrade the critical CRISPR RNA components. Therefore, after neutralizing the trypsin with FBS-containing media, it is crucial to wash the cells with PBS. Alternatively, use “enzyme-free” dissociation media, instead of trypsin.
3. Assemble RNPs individually. For example, use separate reactions for each guide RNA, if targeting multiple sites per sample (e.g., in nickase experiments).
4. We recommend to include the Alt-R Cas9 Electroporation Enhancer in the electroporation. This protocol recommends the use of this non-targeting carrier DNA to improve electroporation efficiency ([Figure 1](#)). Use the same molar concentration of the electroporation enhancer as ribonucleoprotein complex.
5. When using crRNA:tracrRNA duplexes, we recommend using the appropriate Alt-R CRISPR-Cas9 Control Kit for studies in human or mouse cells.

For assistance with control sgRNAs, contact [applicationsupport@idtdna.com](mailto:applicationsupport@idtdna.com).

# CONSUMABLES AND EQUIPMENT

## Consumables from IDT

Kits and reagents	Ordering information
<p><b>Option 1, 2-part guide RNA</b> (crRNA + tracrRNA):</p> <p>Alt-R CRISPR-Cas9 crRNA or Alt-R CRISPR-Cas9 crRNA XT</p> <p>Alt-R CRISPR-Cas9 tracrRNA or Alt-R CRISPR-Cas9 tracrRNA – ATTO 550 or Alt-R CRISPR-Cas9 tracrRNA – ATTO 488 or Alt-R CRISPR-Cas9 tracrRNA – ATTO 647</p>	<p>IDT predesigned and custom crRNA*: <a href="http://www.idtdna.com/CRISPR-Cas9">www.idtdna.com/CRISPR-Cas9</a></p> <p>IDT (cat # 1072532, 1072533, 1072534) IDT (cat # 1075927, 1075928) IDT (cat # 10007810) IDT (cat # 10007853)</p>
<p><b>Option 2, single guide RNA</b> (sgRNA):</p> <p>Alt-R CRISPR-Cas9 sgRNA</p> <p>(Recommended for option 1, 2-part guide RNA)</p> <p>Alt-R CRISPR-Cas9 Control Kit</p> <p>Alt-R <i>S.p.</i> Cas9 Nuclease V3<sup>†</sup></p> <p>Alternatives:</p> <p>Alt-R <i>S.p.</i> HiFi Cas9 Nuclease V3 Alt-R <i>S.p.</i> Cas9 D10A Nickase V3 Alt-R <i>S.p.</i> Cas9 H840A Nickase V3 Alt-R <i>S.p.</i> Cas9 V3, glycerol-free Alt-R <i>S.p.</i> Cas9-GFP V3</p>	<p>IDT predesigned and custom sgRNA*: <a href="http://www.idtdna.com/CRISPR-Cas9">www.idtdna.com/CRISPR-Cas9</a></p> <p>IDT (cat # 1072554 [human] or 1072555 [mouse])</p> <p>IDT (cat # 1081058, 1081059, 10000735)</p> <p>IDT (cat # 1081060, 1081061, 10007803) IDT (cat # 1081062, 1081063) IDT (cat # 1081064, 1081065) IDT (cat# 10007806, 10007807, 10007808) IDT (cat# 10008100, 10008161)</p>
<p>(Optional, but recommended)</p> <p>Alt-R Cas9 Electroporation Enhancer<sup>‡</sup></p>	<p>IDT (cat # 1075915, 1075916, 10007805)</p> <p>Sequence (100 nt): TTAGCTCTGTTTACGTCCCAGCGGGCATGAGAGTAA CAAGAGGGTGTGGTAATATTACGGTACCGAGCACTA TCGATACAATATGTGTCATACGGACACG</p>
<p>Nuclease-Free IDTE, pH 7.5 (1X TE solution)</p>	<p>IDT (cat # 11-01-02-02)</p>

\* We guarantee the performance of our predesigned gRNAs targeting human, mouse, rat, zebrafish, or *C. elegans* genes. For other species, use our proprietary algorithms to design custom gRNAs. If you have protospacer designs of your own or from publications, use our design checker tool to assess their on- and off-targeting potential before ordering gRNAs that are synthesized using our Alt-R gRNA modifications. For details about the predesigned gRNA guarantee, see [www.idtdna.com/CRISPR-Cas9](http://www.idtdna.com/CRISPR-Cas9).

† Alt-R *S.p.* Cas9 Nuclease V3 (wild-type) is suitable for most genome editing studies. However, some experiments may benefit from use of Alt-R *S.p.* HiFi Cas9 Nuclease V3, which has been engineered to reduce off-target effects, while retaining the on-target potency of wild type Cas9. Alt-R Cas9 nickases create single-stranded breaks. When a nickase variant is used with 2 gRNAs, off-target effects are reduced, and homology-directed repair is promoted. Alt-R™ *S.p.* Cas9 V3, glycerol-free may be of interest when working with samples or systems where the presence of glycerol may interfere, such as primary cell cultures or high-throughput instruments with sensitive fluidics.

‡ The enhancer is designed to avoid homology to the human, mouse, or rat genomes, and has been tested as carrier DNA in multiple cell lines, including HEK-293, Jurkat, and K562. Before use in other species, verify that the electroporation enhancer does not have similarity to your host cell genome to limit participation in the double-stranded DNA break repair process.

## Consumables from other suppliers


Kits and reagents	Ordering information
SF Cell Line 96-well Nucleofector™ Kit	Lonza (cat # V4SC-2096)
Dulbecco's Modified Eagle's Medium (DMEM)	ATCC (cat # 30-2002)
Fetal bovine serum (FBS)	General laboratory supplier
Trypsin	General laboratory supplier
1X Phosphate buffered saline (PBS)	General laboratory supplier

## Equipment

Kits and reagents	Ordering information
4D-Nucleofector™ System	Lonza (cat # AAF-1002B with AAF-1002X)
96-well Shuttle™ System	Lonza (cat # AAM-1001S)

# PROTOCOL

## Culture cells

1. Do not use freshly thawed cells for electroporation.
2. Use cells with the lowest passage number possible. Lonza recommends not using HEK-293 cells after passage 20.
3. Replace cell culture media every 2–3 days. For stable cell lines, make sure to include appropriate selection antibiotic.
4. Split cells to maintain confluency  $\leq 90\%$ .
  -  **Note:** Optimal confluency for electroporation of HEK-293 cells with the Nucleofector™ system is 70–85%. Higher cell densities may reduce electroporation efficiency [1].
5. Subculture cells for a minimum of 2–3 days before electroporation, and visually inspect the cells with a microscope to ensure healthy appearance.

## Prepare RNA

1. Resuspend your RNA oligos in IDTE Buffer.

Guide RNA	Final concentration ( $\mu\text{M}$ )
Option 1	
Alt-R CRISPR-Cas9 crRNA	200
Alt-R CRISPR-Cas9 tracrRNA	200
Option 2	
Alt-R CRISPR-Cas9 sgRNA	100

For assistance, use the IDT Resuspension Calculator at [www.idtdna.com/SciTools](http://www.idtdna.com/SciTools).

 **Note:** Store resuspended RNAs at  $-20^{\circ}\text{C}$ .

2. If using sgRNA, proceed to **Form the RNP complex**.
3. Mix the crRNA and tracrRNA oligos in equimolar concentrations in a sterile microcentrifuge tube to a final duplex concentration of  $100\ \mu\text{M}$ . The following table shows an example for a  $10\ \mu\text{L}$  final volume:

Component	Amount ( $\mu\text{L}$ )
$200\ \mu\text{M}$ Alt-R CRISPR-Cas9 crRNA	5
$200\ \mu\text{M}$ Alt-R CRISPR-Cas9 tracrRNA	5
Total volume	10

4. Heat at  $95^{\circ}\text{C}$  for 5 min.
5. Remove from heat and allow to cool to room temperature ( $15\text{--}25^{\circ}\text{C}$ ) on the bench top.

## Form the RNP complex

- For each well undergoing electroporation, dilute the guide RNA and Cas9 enzyme in PBS, gently swirling the pipet tip while pipetting:

Component	Amount ( $\mu\text{L}$ )
PBS	2.1
Alt-R guide RNA (crRNA:tracrRNA duplex from step B5 or sgRNA from step B2)	1.2 (120 pmol)
Alt-R Cas9 enzyme (61 $\mu\text{M}$ stock)*	1.7 (104 pmol)
Total volume	5 <sup>†</sup>

\* All Alt-R *S.p.* Cas9 nucleases and nickases are provided at a stock concentration of 61  $\mu\text{M}$  (10 mg/mL). The GFP-tagged Cas9 is provided at 52  $\mu\text{M}$  (10 mg/mL). Refer to the Application note for tips for using the nickases [4].

<sup>†</sup> The 5  $\mu\text{L}$  reaction volume is for a single nucleofection reaction; scale up as necessary for your experiment.

- Incubate at room temperature for 10–20 min.



**Tip:** To save time, prepare the RNP during the 2 x 10 min centrifugation in **Perform electroporation of cells with Nucleofector™ system**, steps 8 and 11 below.

## Prepare the Nucleofector™ system

- Turn on the Nucleofector™ system and Shuttle device. Open software and log in. Make sure the software connects to the device.
- Open new parameter file.
- Select wells of the 96-well plate, then select the appropriate Nucleofector™ program for your cell line.

## Perform electroporation of cells with the Nucleofector™ system

- Resuspend the Alt-R Cas9 Electroporation Enhancer to 100  $\mu\text{M}$  in IDTE. For assistance, use the Resuspension Calculator at [www.idtdna.com/SciTools](http://www.idtdna.com/SciTools).
- Add the entire Supplement to the Nucleofection™ Solution SF before first use, as directed by the manufacturer.
- Prepare a 96-well culture plate to receive cells following nucleofection.




**Note:** We recommend dividing each nucleofection into 3 replicate wells.

- Fill necessary wells with 175  $\mu\text{L}$  of culture media (DMEM, 10% FBS).
  - Store in a tissue culture incubator (37°C, 5% CO<sub>2</sub>).
- Prepare an additional aliquot of 75  $\mu\text{L}$ /well of culture media (DMEM, 10% FBS) and pre-warm to 37°C.
  - Harvest cells in a 150 cm<sup>2</sup> flask by trypsinization.
    - Aspirate media from cultured cells and wash once with 10 mL of 1X PBS.
    - Add 4 mL of 1X trypsin solution to cells and incubate at 37°C until cells just release from flask.
 

**Important:** Avoid over-trypsinizing cells while making sure that you achieve a single-cell suspension.
    - Neutralize trypsin by adding 16 mL of culture media (DMEM, 10% FBS).
  - Count the trypsinized cells.



7. Transfer the total number of cells you need for your experiment to a sterile 15 mL tube.

 **Note:** Typically,  $2 \times 10^5$ – $5 \times 10^5$  cells per well is the optimal range of HEK-293 cells for nucleofection [1]. In developing this protocol we used  $3.5 \times 10^5$  HEK-293 cells per nucleofection; scale up for the appropriate number of wells.

8. Centrifuge the cells at 30 x g for 10 min at room temperature.

9. Remove as much supernatant as possible without disturbing the pellet.

 **Important!** Wash cells in 5 mL of 1X PBS.

Trypsin and FBS may contain RNase activity. Therefore, after neutralizing trypsin with FBS-containing media, we recommend washing the cells with PBS. A good alternative would be to use “enzyme-free” dissociation media.

10. Centrifuge at 30 x g for 10 min at room temperature.


11. Remove as much supernatant as possible without disturbing the pellet.

12. Resuspend cells by adding 20  $\mu$ L of supplemented Nucleofector™ Solution SF (from [Perform electroporation of cells with Nucleofector™ system](#), step 1) per  $3.5 \times 10^5$  cells.

13. Pipet 20  $\mu$ L of cell suspension into each well of a V-bottom plate.

14. To each well, add 5  $\mu$ L of the RNP (from [Form the RNP complex](#), step 2) and 1  $\mu$ L of 100  $\mu$ M Alt-R Cas9 Electroporation Enhancer (from [Perform electroporation of cells with Nucleofector™ system](#), step 1).


15. Pipet to mix, then transfer 25  $\mu$ L of the cell:RNP complex mixture to the wells of the 96-well Nucleocuvette™ module.

 **Note:** The total volume is 26  $\mu$ L, containing 4.6  $\mu$ M guide RNA, 4  $\mu$ M Cas9 nuclease, and 4  $\mu$ M Cas9 electroporation enhancer. The total volume allows for easy pipetting of the 25  $\mu$ L required volume in the electroporation cuvette, and prevents the formation of air bubbles.

16. Gently tap the Nucleocuvette™ module to make sure no air bubbles are present.

17. Place Nucleocuvette™ module in Shuttle device, and select **Upload and start**.

18. After electroporation, remove the Nucleocuvette™ module from the instrument.

 **Note:** Be sure to save the file before nucleofection occurs.

19. Add 75  $\mu$ L of pre-warmed culture media (from [Perform electroporation of cells with Nucleofector™ system](#), step 4) per well and resuspend cells by gently pipetting up and down.

20. Transfer 25  $\mu$ L of resuspended cells to the 175  $\mu$ L of culture media (DMEM, 10% FBS) from [Perform electroporation of cells with Nucleofector™ system](#), step 3 in triplicate.

21. Incubate cells in a tissue culture incubator (37°C, 5% CO<sub>2</sub>) for 48 hr.

To investigate on-target mutations with the mismatch endonuclease T7EI, use the [Alt-R Genome Editing Detection Kit](#) (cat # 1075931, 1075932, 1075933) [5].

## REFERENCES

---

1. Lonza. (2010) **Amaxa™ 4D-Nucleofector™ Optimization Protocol for Cell Lines** [Online] Basel, Switzerland, Lonza, Ltd. [Accessed 27 Aug 2021]
2. Lonza. (2009) **Amaxa 96-well Shuttle Protocol for HEK-293 (ATCC®)**. [Online] Basel, Switzerland, Lonza, Ltd. [Accessed 6 Jun 2021]
3. Integrated DNA Technologies. (2017) **Alt-R CRISPR-Cas9 System: User Guide for cationic lipid delivery of CRISPR-Cas9 ribonucleoprotein into mammalian cells**. [Online] Coralville, IA, Integrated DNA Technologies, Inc. [Accessed 6 Jun 2021]
4. Yan S, Schubert M, et al. (2017) **Applications of Cas9 nickases for genome engineering**. [Online] Coralville, IA, Integrated DNA Technologies, Inc. [Accessed 6 Jun 2021]
5. Integrated DNA Technologies. (2017) **Alt-R Genome Editing Detection Kit protocol**. [Online] Coralville, IA, Integrated DNA Technologies, Inc. [Accessed 6 Jun 2021]

## Alt-R CRISPR-Cas9 System: Delivery of ribonucleoprotein complexes into HEK-293 cells using the Lonza® Nucleofector™ System

Technical support: [applicationsupport@idtdna.com](mailto:applicationsupport@idtdna.com)

For more than 30 years, IDT's innovative tools and solutions for genomics applications have been driving advances that inspire scientists to dream big and achieve their next breakthroughs. IDT develops, manufactures, and markets nucleic acid products that support the life sciences industry in the areas of academic and commercial research, agriculture, medical diagnostics, and pharmaceutical development. We have a global reach with personalized customer service.

> SEE WHAT MORE WE CAN DO FOR YOU AT [WWW.IDTDNA.COM](http://WWW.IDTDNA.COM).

**For Research Use Only. Not intended for use in diagnostic procedures.** Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2021 Integrated DNA Technologies, Inc. All rights reserved. Alt-R and Ultramer are trademarks of Integrated DNA Technologies, Inc., and registered in the USA. ATTO is a trademark of ATTO-TEC GmbH. All other marks are the property of their respective owners. For specific trademark and licensing information, see [www.idtdna.com/trademarks](http://www.idtdna.com/trademarks).  
MAPPS Doc ID# RUO21-0212\_001 09/21