

For resuspending and annealing dried-down, complementary oligonucleotides to form a duplex

1. Centrifuge tubes before opening to ensure duplexed oligos are at the bottom of the tube.
2. Resuspend both oligos in Nuclease-Free Duplex Buffer* (Cat # 11-01-03-01) to reach the appropriate final volume. For example:

	100 μ M final duplex concentration		50 μ M final duplex concentration		10 μ M final duplex concentration	
Oligo amount†	Buffer volume for oligos	Total volume for annealing	Buffer volume for oligos	Total volume for annealing	Buffer volume for oligos	Total volume for annealing
10 nmol	50 μ L each	100 μ L	100 μ L each	200 μ L	500 μ L each	1 mL
25 nmol	125 μ L each	250 μ L	250 μ L each	500 μ L	1.5 mL each	2.5 mL
50 nmol	250 μ L each	500 μ L	500 μ L each	1 mL	2.5 mL each	5 mL

† Refer to the IDT resuspension calculator at www.idtdna.com to calculate dilutions for other nanomole amounts.

3. Mix the two oligos in equimolar concentrations.
4. Heat at 94°C for 2 min.
5. Remove from heat and allow to cool to room temperature.
6. If needed, dilute the annealed oligonucleotides using Nuclease-Free Duplex Buffer or 1X IDTE Buffer (Cat # 11-01-02-05).
7. Store at -20°C.

 **Note:** If the annealed oligos are to be used on multiple occasions, divide into smaller aliquots and store at -20°C.

* Nuclease-Free Duplex Buffer is certified nuclease-free by testing with RNaseAlert® and DNaseAlert™ reagents from IDT. For more information visit www.idtdna.com.

For Research Use Only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2022 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks. Doc ID: RUO22-1160_001 9/22



* 1 0 7 8 2 5 2 *