

## TA cloning method

For use with IDT Gene Fragments, such as:

- gBlocks™ Gene Fragments
- gBlocks HiFi Gene Fragments
- eBlocks™ Gene Fragments

> SEE WHAT MORE WE CAN DO FOR YOU AT [WWW.IDTDNA.COM](http://WWW.IDTDNA.COM).

NGS • CRISPR • Functional Genomics • Synthetic Biology • PCR & qPCR • Custom Oligos • Contract Manufacturing

## Table of contents

Introduction	3
Guidelines	4
Advantages of the TA cloning method	4
Limitations to the TA cloning method	4
IDT Gene Fragments	4
gBlocks Gene Fragments	4
gBlocks HiFi Gene Fragments	4
eBlocks Gene Fragments	4
Consumables and equipment	5
Consumables—IDT	5
Consumables—Other suppliers	5
Equipment	6
Protocol	7
Resuspend your IDT Gene Fragment	7
Store your Gene Fragments	7
gBlocks & gBlocks HiFi Gene Fragments	7
eBlocks Gene Fragments	7
A-tail your Gene Fragment	8
Prepare vector	8
Option 1: Linearize by restriction digestion	8
Option 2: Linearize by amplification and digest with DpnI	9
Add T-tail to vector with terminal transferase	10
Clone gene fragment into vector	10
Transform	11
Example protocol	11
References	11


# INTRODUCTION

TA cloning is a fast and relatively simple way to clone relatively short pieces of DNA with moderate efficiency. Several commercially available cloning kits use the TA cloning method. Vectors in these kits have single "T" base overhangs at the vector cloning site that is designed to complement the residual 3' "A" base left by Taq polymerase during PCR reactions. To TA clone DNA with blunt ends using a product such as IDT Gene Fragments, "A" base overhangs will be added after a brief incubation with Taq polymerase in the presence of dATP (**Figure 1**).

## TA cloning

"A" tail insert with Taq polymerase

Ligate insert into vector



**Figure 1. Overview of TA cloning.**

Compared to methods such as seamless cloning, TA cloning is non-directional and as such, the method cannot be used to assemble more than one fragment at a time. Nevertheless, it can be a good cloning option to rapidly clone PCR products for routine sequence identification, subcloning applications, and creation of small recombinant plasmids.

This protocol includes how to prepare your blunt-ended gene fragment for TA cloning. If you are using a TA cloning kit, follow the manufacturer's instructions for cloning your prepared fragment. If you are not using a TA cloning kit, this protocol will provide guidelines for preparing your TA cloning vector, ligating your fragment into the vector, and transforming your plasmid into bacteria for selection and screening.

# GUIDELINES

---

## Advantages of the TA cloning method

- More efficient than blunt-end cloning
- Restriction enzymes are not required to ligate the insert with the vector

## Limitations to the TA cloning method

- Not directional—inserts will be cloned in both the forward and the reverse orientation
- Inserts an additional “A” and “T” base at the cloning junctions
- Sample must be purified to avoid selective cloning of truncated products
- Blunt inserts not generated by traditional PCR methods require A-tailing before cloning
- Requires specially prepared plasmids with “T” base overhangs
- Inefficient when cloning more than one fragment

## IDT Gene Fragments

### gBlocks Gene Fragments

gBlocks Gene Fragments are synthetic, double-stranded DNA (125–3000 bp), delivered dry and normalized to 250, 500, or 1000 ng, depending on length. These fragments are ideal for gene assembly, genome editing, qPCR controls, and more.

### gBlocks HiFi Gene Fragments

gBlocks HiFi Gene Fragments are between 1000–3000 bp in length, shipped dry in tubes, and normalized at 1000 ng. These high-fidelity fragments are intended for the assembly of large constructs.

### eBlocks Gene Fragments

eBlocks Gene Fragments are synthetic, double-stranded DNA (300–1500 bp), normalized to 200 ng and delivered at 10 ng/μL in Nuclease-Free Water. These gene fragments are uniquely suited for high-throughput screening and can ship in as little as 1–3 business days from order confirmation.

# CONSUMABLES AND EQUIPMENT

## Consumables—IDT

Item	Catalog #
Gene fragments, such as: gBlocks Gene Fragments gBlocks HiFi Gene Fragments eBlocks Gene Fragments	<a href="http://www.idtdna.com/gBlocks">www.idtdna.com/gBlocks</a> <a href="http://www.idtdna.com/gBlocks">www.idtdna.com/gBlocks</a> <a href="http://www.idtdna.com/eBlocks">www.idtdna.com/eBlocks</a>
Nuclease-Free Water	11-04-02-01

## Consumables—Other suppliers

Item
TE buffer
Taq DNA polymerase
1 mM dATP
25 mM MgCl <sub>2</sub>
TA cloning vector kit*
If not using a TA cloning vector kit: (Option 1) To linearize vector manually by restriction digestion: <ul style="list-style-type: none"><li>• Vector</li><li>• Restriction enzyme (blunt cutting enzyme, such as EcoRV)</li><li>• Alkaline phosphatase</li></ul> (Option 2) To linearize vector manually by PCR amplification: <ul style="list-style-type: none"><li>• Vector</li><li>• Forward and reverse primers</li><li>• dNTPs</li><li>• High fidelity polymerase</li><li>• MgSO<sub>4</sub></li><li>• DpnI</li></ul>
To A-tail and ligate linearized vector: <ul style="list-style-type: none"><li>• 2.5 mM CaCl<sub>2</sub></li><li>• 1 mM ddTTP</li><li>• Terminal transferase</li><li>• T4 ligase</li></ul>
Cell transformation reagents <ul style="list-style-type: none"><li>• Competent bacteria</li><li>• 14 mL BD Falcon® polypropylene tube (Corning)</li><li>• SOC media</li><li>• LB plates</li></ul>

\* There are several commercial TA cloning solutions available from various manufacturers; follow the protocols for the TA cloning vector or kit to ensure cloning success.

## Equipment

Item
PCR thermal cycler
(Optional) Heat blocks at 37°C and 80°C
Water bath at 42°C
Incubator at 37°C

# PROTOCOL

## Resuspend your IDT Gene Fragment

**!** **Important:** The dried down IDT Gene Fragment pellet can become displaced from the bottom of the tube during shipping, so it is important to centrifuge your tubes or plates before opening them.

1. Centrifuge the tube for 3–5 seconds at a minimum of 3000 x g to pellet the material to the bottom of the tube.
2. Add TE buffer to the tube to reach your desired final concentration.

Final stock concentration	Resuspension volume of TE buffer (µL) for IDT Gene Fragments		
	Synthesis scale		
	250 ng	500 ng	1000 ng
10 ng/µL	25	50	100
20 ng/µL	Not recommended	25	50
50 ng/µL	Not recommended	10	20

**☰** **Note:** Gene fragments can be resuspended in Nuclease-Free Water in some cases (see [Store your Gene Fragments](#)).

3. Briefly vortex, then centrifuge.

## Store your Gene Fragments

### gBlocks & gBlocks HiFi Gene Fragments

Store gBlocks and gBlocks HiFi Gene Fragments at –20°C in TE up to 24 months. For short-term storage of less than one month, they can be stored at –20°C in Nuclease-Free Water. Aliquot resuspended gene fragments to avoid more than 2–3 freeze-thaw cycles.

### eBlocks Gene Fragments

eBlocks fragments are delivered at 10 ng/µL in Nuclease-Free Water. Store eBlocks Gene Fragments at –20°C up to 24 months. For short-term storage of less than two weeks, store at 4°C.

## A-tail your Gene Fragment

The tailing procedure modifies blunt-ended IDT Gene Fragments with a single “A” base overhang on the 3' ends, resulting in DNA that is compatible with TA cloning kits.

1. Add these reagents together for a 20  $\mu$ L reaction:

Reagent	Amount
IDT Gene Fragment	50 ng
Taq polymerase	1–3 units
10X Taq polymerase buffer minus Mg	2 $\mu$ L
1 mM dATP	To [0.2 mM]
25 mM MgCl <sub>2</sub>	To [1.5 mM]
Nuclease-Free Water	To final 20 $\mu$ L volume
<b>Total volume</b>	<b>20 <math>\mu</math>L</b>


2. Incubate at 70°C for 15–30 minutes.
3. Depending on if you are using a kit or not, follow one of these two steps:
  - When using a kit, use 1–10  $\mu$ L of the A-tailed gene fragment for TA cloning according to the manufacturer's instructions.
  - When not using a kit, proceed further into the protocol for guidance on cloning your insert into your vector.

## Prepare vector

Vectors used for TA cloning can be purchased commercially, or as described here, prepared manually in a two-step process [1]. The first step is to linearize the vector with a blunt cutting restriction endonuclease digestion (**Option 1**), or with PCR (**Option 2**). The second step is to tail the linear vector with dideoxythymidine triphosphate (ddTTP) using terminal transferase.

### Option 1: Linearize by restriction digestion

Supercoiled vector isolated from *E. coli* or purchased from a commercial vendor can be linearized using a blunt cutting restriction enzyme, provided a restriction site is present in the vector. In addition, the linearized vector should then be dephosphorylated using a phosphatase to prevent religation of the empty vector's ends. This protocol provides an example using EcoRV and Thermosensitive Alkaline Phosphatase.

 **Note:** The volumes and concentrations are examples; follow the manufacturer's instructions for materials specific to your application.


Add these reagents together, then incubate at 37°C for 1 hour, followed by 80°C for 20 minutes.

Reagent	Amount
Plasmid	10 $\mu$ g
10X Buffer #3	2.5 $\mu$ L
EcoRV (400 U/ $\mu$ L)	1 $\mu$ L
BSA	0.5 $\mu$ L
Nuclease-Free Water	To final 25 $\mu$ L volume
<b>Total volume</b>	<b>25 <math>\mu</math>L</b>



## Option 2: Linearize by amplification and digest with DpnI

Alternatively, vectors can be amplified using primers that have their 5' ends at the insertion site (i.e., oriented to extend away from the insertion site). For amplification, use a high-fidelity polymerase that leaves blunt ends on the products. To remove the PCR template, digest the reaction using DpnI, which will only digest Dam-methylated DNA isolated from *E. coli* and not the PCR-amplified DNA.

 **Note:** The volumes, temperatures, and concentrations are only examples; follow the manufacturer's instructions for materials specific to your application.


1. Set up the components (**Table 1**) for each of your vectors.

Table 1. Components of PCR.

Reagent	Amount
Supercoiled plasmid	1 ng
5 $\mu$ M Forward Primer	1 $\mu$ L
5 $\mu$ M Reverse Primer	1 $\mu$ L
2 mM dNTPs	2.5 $\mu$ L
10X KOD Buffer	2.5 $\mu$ L
25 mM MgSO <sub>4</sub>	1.5 $\mu$ L
High-fidelity polymerase (2.5 U/ $\mu$ L)	0.5 $\mu$ L
Nuclease-Free Water	To final 25 $\mu$ L volume
<b>Total volume</b>	<b>25 <math>\mu</math>L</b>

2. (Optional) Confirm the linear product was generated by running 5  $\mu$ L on a 0.8% agarose gel with a DNA ladder and 200 ng of uncut plasmid.
3. Digest the template from the PCR-amplified vector with DpnI by incubating the following reagents at 37°C for 1 hour:

Reagent	Amount
PCR product	17 $\mu$ L
10X Buffer	2 $\mu$ L
DpnI, 20 U/ $\mu$ L	1 $\mu$ L
<b>Total volume</b>	<b>20 <math>\mu</math>L</b>

 **Tip:** Standard PCR primers do not contain 5' phosphates, so you do not need to dephosphorylate after PCR.

4. (Optional) Purify the PCR-amplified vector using a PCR cleanup kit.
5. Confirm the linear product was generated by running 5  $\mu$ L on a 0.8% agarose gel with a DNA ladder and 200 ng of uncut plasmid.

## Add T-tail to vector with terminal transferase

1. Create the following reaction setup, noting that the molar amount of plasmid can be calculated with the following equation:

Mass of vector that equals 5 pmol (in  $\mu\text{g}$ ) = length of plasmid in bp \* 650 Daltons \*  $5 \times 10^{-6}$

 **Note:** A 2700 bp plasmid like pUC19 requires approximately 8.8  $\mu\text{g}$  of material.

Reagent	Amount
10X terminal transferase buffer	5 $\mu\text{L}$
2.5 mM $\text{CaCl}_2$ solution	5 $\mu\text{L}$
Linear plasmid DNA	5 pmol
1 mM ddTTP	0.5 $\mu\text{L}$
Terminal transferase	10 units
Nuclease-Free Water	To final 40 $\mu\text{L}$ volume
<b>Total volume</b>	<b>40 <math>\mu\text{L}</math></b>

2. Incubate the reaction for 20 minutes at 37°C.
3. Stop the reaction by heating to 70°C for 10 minutes.

## Clone gene fragment into vector

T4 ligases are available in standard or “quick” versions. In general, the quick versions contain a crowding agent such as PEG that increases the ligation kinetics and decreases reaction time. Quick versions are preferred for day-to-day uses.

To ligate your IDT Gene Fragment into the vector efficiently, use the optimum molar ratio of vector to gene fragment. The ratio we recommend is 1:5–1:12 vector to gene fragment.

### Tips:

- High salt concentrations from DNA can inhibit ligation activity, so the use of DNA cleanup kits can often improve reaction efficiencies.
  - T4 DNA ligases require buffers containing ribo ATP. Always use fresh buffers that have not been repeatedly frozen and thawed to ensure the ATP is active.
  - Be careful to not heat inactivate ligase reactions with buffers containing PEG as this decreases transformational efficiency.
4. Set up the reaction by combining these items:

Item	Amount
Linearized vector	50 ng
IDT Gene Fragment	5–12X molar excess*
2X ligase buffer	10 $\mu\text{L}$
T4 DNA ligase	1 $\mu\text{L}$
Nuclease-Free Water	Up to 20 $\mu\text{L}$
<b>Total volume</b>	<b>20 <math>\mu\text{L}</math></b>

\* Molar ratios of the gBlocks, gBlocks HiFi, or eBlocks Gene Fragments (IDT) can be converted to ng using the following formula:

$$50 \text{ ng} \times \text{desired molar ratio} \times \frac{\text{IDT Gene Fragments length (bp)}}{\text{Plasmid length (bp)}} = \text{ng gene fragment needed}$$

5. (Optional) Store at –20°C.

## Transform

Several lines of competent *E. coli* can be purchased from a variety of vendors and provide a reliable way to achieve high transformation efficiencies. Alternatively, competent cells can be prepared in the lab by following the protocols outlined in Sambrook et al. [2].

### Example protocol

Here is an example protocol for transformation; follow the manufacturer's instructions for materials specific to your product.

1. Thaw cells on wet ice.
2. Add 25  $\mu$ L cells to a pre-chilled 14 mL BD Falcon polypropylene tube on ice.
3. Add 2  $\mu$ L of ligation mixture and mix gently.
4. Incubate on wet ice for 30 minutes.
5. Place in a 42°C water bath for 45 seconds.
6. Return to ice for 2 minutes.
7. Add 250  $\mu$ L of SOC media to the cells, then incubate shaking at 37°C for 1 hour.
8. Plate 125  $\mu$ L on LB plates with the appropriate selection reagents for your vector.
9. Incubate the plates inverted in a 37°C incubator overnight.
10. Select and screen several colonies.

The resulting completed plasmid is ready for sequencing.

## REFERENCES

---

1. Holton TA, Graham MW. **A simple and efficient method for direct cloning of PCR products using ddT-tailed vectors**. *Nucleic Acids Res.* 1991; 19(5):1156. doi:10.1093/nar/19.5.1156
2. Green M, Sambrook J. **Molecular Cloning: A Laboratory Manual**. (Fourth Edition). Cold Spring Harbor Laboratory Press; 2012.

## TA cloning method

Contact Us: [www.idtdna.com/contactus](http://www.idtdna.com/contactus)

For more than 30 years, IDT's innovative tools and solutions for genomics applications have been driving advances that inspire scientists to dream big and achieve their next breakthroughs. IDT develops, manufactures, and markets nucleic acid products that support the life sciences industry in the areas of academic and commercial research, agriculture, medical diagnostics, and pharmaceutical development. We have a global reach with personalized customer service.

> SEE WHAT MORE WE CAN DO FOR YOU AT [WWW.IDTDNA.COM](http://WWW.IDTDNA.COM).

**For Research Use Only. Not for use in diagnostic procedures.**

Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2022 Integrated DNA Technologies, Inc. All rights reserved. gBlocks is a trademark of Integrated DNA Technologies, Inc. and is registered in the USA. eBlocks is a trademark of Integrated DNA Technologies, Inc. All other marks are the property of their respective owners. For specific trademark and licensing information, see [www.idtdna.com/trademarks](http://www.idtdna.com/trademarks).  
Doc ID: RUO22-1132\_001 11/22