

Follow the instructions and incubations for each step. All mixing steps should be performed on ice. Pipette up and down 8 times or vortex to mix after re-suspending each lysphere and spin down prior to incubations and transfers. For incubations, use a lid heated  $\geq 100^{\circ}\text{C}$ , except where specified otherwise.

<b>Random Priming</b>	Add 20-250 ng RNA to Random Priming reagent for a total volume of 20 $\mu\text{L}$ .  <b>Incubate as indicated</b> 🖱️	<table border="1"> <thead> <tr> <th>Step</th> <th>Incubation Temperature</th> <th>Incubation Time</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>65°C</td> <td>5 min</td> </tr> <tr> <td>2</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table>	Step	Incubation Temperature	Incubation Time	1	65°C	5 min	2	4°C	Hold						
	Step	Incubation Temperature	Incubation Time														
1	65°C	5 min															
2	4°C	Hold															
<b>First Strand cDNA Synthesis</b>	Transfer 20 $\mu\text{L}$ Random Priming mixture to First Strand cDNA Synthesis reagent.  <b>Incubate as indicated</b> 🖱️	<table border="1"> <thead> <tr> <th>Step</th> <th>Incubation Temperature</th> <th>Incubation Time</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>25°C</td> <td>10 min</td> </tr> <tr> <td>2</td> <td>42°C</td> <td>30 min</td> </tr> <tr> <td>3</td> <td>80°C</td> <td>20 min</td> </tr> <tr> <td>4</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table>	Step	Incubation Temperature	Incubation Time	1	25°C	10 min	2	42°C	30 min	3	80°C	20 min	4	4°C	Hold
	Step	Incubation Temperature	Incubation Time														
	1	25°C	10 min														
	2	42°C	30 min														
3	80°C	20 min															
4	4°C	Hold															
Transfer 1 $\mu\text{L}$ First Strand cDNA Synthesis mixture to 9 $\mu\text{L}$ water for PreSeq RNA QC Assay.																	
<b>Second Strand cDNA Synthesis</b>	Add 21 $\mu\text{L}$ ultrapure water to 19 $\mu\text{L}$ First Strand cDNA Synthesis mixture and then transfer to Second Strand cDNA Synthesis reagent.  <b>Incubate as indicated</b> 🖱️	<table border="1"> <thead> <tr> <th>Step</th> <th>Incubation Temperature</th> <th>Incubation Time</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>16°C</td> <td>60 min</td> </tr> <tr> <td>2</td> <td>75°C</td> <td>20 min</td> </tr> <tr> <td>3</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table>	Step	Incubation Temperature	Incubation Time	1	16°C	60 min	2	75°C	20 min	3	4°C	Hold			
	Step	Incubation Temperature	Incubation Time														
	1	16°C	60 min														
	2	75°C	20 min														
3	4°C	Hold															
<b>Optional stopping point after this step. Store at -10°C to -30°C.</b>																	
<b>PreSeq RNA QC Assay</b>	Build qPCR reaction and incubate as per your master mix-specific instructions.  iTaq SYBR Green Supermix (Bio-Rad) – 5 $\mu\text{L}$ 10X VCP Primer Mix – 1 $\mu\text{L}$ Diluted cDNA sample or NTC – 4 $\mu\text{L}$  <b>Incubate as indicated</b> 🖱️	<table border="1"> <thead> <tr> <th>Incubation Temperature</th> <th>Incubation Time</th> <th># of cycles</th> </tr> </thead> <tbody> <tr> <td>95°C</td> <td>20 [20#] sec</td> <td>1</td> </tr> <tr> <td>95°C</td> <td>3 [15#] sec</td> <td rowspan="2">35</td> </tr> <tr> <td>60°C</td> <td>30 [60#] sec</td> </tr> <tr> <td>60-95°C</td> <td>0.5°C/sec increment</td> <td>1</td> </tr> </tbody> </table> <p># Times in [] are for standard cycling</p>	Incubation Temperature	Incubation Time	# of cycles	95°C	20 [20#] sec	1	95°C	3 [15#] sec	35	60°C	30 [60#] sec	60-95°C	0.5°C/sec increment	1	
	Incubation Temperature	Incubation Time	# of cycles														
	95°C	20 [20#] sec	1														
	95°C	3 [15#] sec	35														
60°C	30 [60#] sec																
60-95°C	0.5°C/sec increment	1															
<b>End Repair</b>	Transfer 40 $\mu\text{L}$ Second Strand cDNA Synthesis mixture to End Repair reagent.  <b>Incubate as indicated with <u>unheated lid</u></b> 🖱️	<table border="1"> <thead> <tr> <th>Step</th> <th>Incubation Temperature</th> <th>Incubation Time</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>25°C</td> <td>30 min</td> </tr> <tr> <td>2</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table>	Step	Incubation Temperature	Incubation Time	1	25°C	30 min	2	4°C	Hold						
	Step	Incubation Temperature	Incubation Time														
	1	25°C	30 min														
2	4°C	Hold															
	AMPure™ XP clean-up (100 $\mu\text{L}$ ). Elute in 20 $\mu\text{L}$																

**Ligation Step 1**

Transfer 20 µL End Repair mixture to Ligation Step 1 reagent.

**Incubate as indicated** 🖱️

Step	Incubation Temperature	Incubation Time
1	37°C	15 min
2	4°C	Hold

AMPure™ XP clean-up (50 µL). Elute in 42 µL

**MBC Adapters**

Transfer 40 µL Ligation Step 1 mixture to the MBC adapters.

**Ligation Step 2**

Transfer the entire volume of the MBC adapters to Ligation Step 2 reagent.

**Incubate as indicated with unheated lid** 🖱️



**Optional stopping point before purification. Store at -10°C to -30°C.**

Step	Incubation Temperature	Incubation Time
1	22°C	5 min
2	4°C	Hold

Ligation clean-up beads (50 µL). Elute in 18 µL 5 mM NaOH, 75°C 10 min

**First PCR**

Add 2 µL GSP1 primers to the First PCR reagent. Then add 18 µL of purified Adapter Ligation mixture and mix.

**Incubate as indicated** 🖱️



**Optional stopping point after this step. Store at -10°C to -30°C.**

Incubation Temperature	Incubation Time	# of cycles
95°C	3 min	1
95°C	30 sec	15*
65°C *	5 min* (ramp rate 100%)	
72°C	3 min	1
4°C	Hold	1

AMPure™ XP clean-up (24 µL). Elute in 20 µL

**Second PCR**

Add 2 µL GSP2 primers to Second PCR reagent. Then add 18 µL of purified First PCR mixture and mix.

**Incubate as indicated** 🖱️



**Optional stopping point after this step. Store at -10°C to -30°C.**

Incubation Temperature	Incubation Time	# of cycles
95°C	3 min	1
95°C	30 sec	20*
65°C *	5 min* (ramp rate 100%)	
72°C	3 min	1
4°C	Hold	1

AMPure™ XP clean-up (24 µL). Elute in 20 µL

\* Refer to Product Insert for panel-specific parameters

**Proceed with Protocol: Quantify, Normalize, and Sequence Protocol for Illumina®**