


Pooling multiple rhAmpSeq panels

When pooling multiple rhAmpSeq panels, first determine which of the following scenarios applies to you:

- 1 All panels include individual oligos in plates.
- 2 All panels include preformulated oligos in tubes.
- 3 Panels include both preformulated tubes and individual oligos in plates.

Figure 1. Overview of scenarios when pooling multiple rhAmpSeq panels.



* You must redesign your panel using the 'Make compatible with a pre-existing panel' option in the [rhAmpSeq Design Tool](#).

Scenario 1: Combining panels with individual oligos in plates

1. Resuspend each oligo (rhAmp primer) to 50 μM , then make separate forward and reverse pools by combining equal volumes of each primer from all panels.

- a. Use the following formula to determine the amount of IDTE, pH 7.5 needed to resuspend each to a final concentration of 50 μM (0.05 nmol/ μL):

$$\text{Amount of IDTE, pH 7.5} = \frac{\text{nmol of rhAmp primer}}{0.05 \text{ nmol}/\mu\text{L}}$$

For example:

$$\text{Amount of IDTE, pH 7.5} = \frac{0.4 \text{ nmol of rhAmp primer}}{0.05 \text{ nmol}/\mu\text{L}} = 8 \mu\text{L}$$

! **Important!** Use IDTE, pH 7.5 (10 mM Tris; 0.1 mM EDTA, pH 7.5), or equivalent. Do not resuspend in water or TE buffer.

- b. Centrifuge all rhAmp primer plates before resuspension.
- c. Resuspend each rhAmp primer to 50 μM using the volume of IDTE, pH 7.5 determined in step a.
- d. Seal the plates, vortex to fully resuspend, then briefly centrifuge.
2. Combine equal volumes of each 50 μM stock rhAmp primer to create stock forward or reverse primer pools.

☰ **Note:** Forward primers are combined to make a forward rhAmp Primer Pool, and reverse primers are combined to make a reverse rhAmp Primer Pool.

! **Important!** Do not combine forward and reverse rhAmp primers until you are ready to perform Targeted rhAmpSeq PCR 1 in the rhAmpSeq library preparation standard protocol [1], or the rhAmpSeq high-throughput library preparation protocol [2].

3. Calculate the total number of primers in each pool.
 - a. If there are <500 primers, refer to [Table 1](#) and dilute each pool to the appropriate concentration using IDTE, pH 7.5 to create a 10X concentration.

Table 1. 10X primer pool calculations.

Panel size	10X rhAmp primer concentration
≥500-plex	50 μM (total)
101-plex through 499-plex	100 nM (each primer)
≤100-plex	250 nM (each primer)

- b. If there are ≥500 primers, no dilution is needed because your pools are already at a 10X concentration.
4. Proceed directly to “Perform Targeted rhAmpSeq PCR 1” in your protocol [1] or [2].

Example 1

Panel 1: 250-plex, forward (FWD) and reverse (REV) pools ordered at 0.4 nmol scale in plates

Panel 2: 20-plex, FWD and REV pools ordered at 0.4 nmol scale in plates

Final pooled panel size: 270 primers in each forward and reverse pool

Final 10X rhAmp primer concentration (see Table 1): 100 nM each primer

- Resuspend dried down oligos to 50 μM:
Add 8 μL of IDTE, pH 7.5 buffer to each well.
- Create a single forward pool and separate reverse pool:
 - Combine 2 μL of forward primers (from Panel 1 and Panel 2) to create the forward pool.
 - Combine 2 μL of reverse primers (from Panel 1 and Panel 2) to create the reverse pool.

Each pool has a total volume of 540 μL.



Note: If the pooled panel had been ≥500-plex, no further dilution of the panel would be needed; therefore, skip to step e.

- Calculate individual primer concentration:
 $50 \mu\text{M} / 270 \text{ primers} = 185 \text{ nM per primer}$.
- Dilute each pool to 100 nM per primer:
Combine 100 μL of stock pool (185 nM) with 85 μL of IDTE, pH 7.5.

The total volume of each forward and reverse pool should be 185 μL (100 nM per primer).
- Proceed directly to “Perform Targeted rhAmpSeq PCR 1” in your protocol [1] or [2].

Example 2

Panel 1: 80-plex, FWD and REV pools ordered at 4 nmol scale in plates

Panel 2: 50-plex, FWD and REV pools ordered at 0.4 nmol scale in plates

Final pooled panel size: 130 primers in each forward and reverse pool

Final 10X rhAmp primer concentration (see Table 1): 100 nM each primer

- Resuspend dried down oligos to 50 μM:
 - Panel 1: Add 80 μL of IDTE, pH 7.5 buffer to each well.
 - Panel 2: Add 8 μL of IDTE, pH 7.5 buffer to each well.
- Create a single forward pool and separate reverse pool:
 - Combine 2 μL of forward primers (from Panel 1 and Panel 2) to create the forward pool.
 - Combine 2 μL of reverse primers (from Panel 1 and Panel 2) to create the reverse pool.

Each pool has a total volume of 260 μL.



Note: If the pooled panel had been ≥500-plex, no further dilution of the panel would be needed; therefore, skip to step e.

- Calculate individual primer concentration:
 $50 \mu\text{M} / 130 \text{ primers} = 385 \text{ nM per primer}$.
- Dilute each pool to 100 nM per primer:
Combine 100 μL of stock pool (385 nM) with 285 μL of IDTE, pH 7.5.

The total volume of each forward and reverse pool should be 385 μL (100 nM per primer).
- Proceed directly to “Perform Targeted rhAmpSeq PCR 1” in your protocol [1] or [2].

Scenario 2: Combining preformulated oligos in tubes

1. Calculate the individual primer concentration in each panel.
 - a. ≥ 500 -plex panels: Divide 50 μM by the number of primers in the pool (i.e., 50 $\mu\text{M}/1000$ primers = 50 nM/primer).
 - b. ≤ 499 -plex panels: Use [Table 1](#) to determine that individual primer concentration in each tube.
2. Determine which of the panels has the lower per primer concentration, then dilute other panels to that concentration.
3. Combine equal volumes of each panel to create separate forward and reverse pooled panels, then proceed directly to "Perform Targeted rhAmpSeq PCR 1" in your protocol [1] or [2].



Note: Resulting per primer concentration will be lower than amount recommended in [Table 1](#). Reaction and thermal cycling parameters may require optimization for optimal assay performance.

Example 1

Panel 1: 250-plex, forward (FWD) and reverse (REV) pools provided as preformulated oligos in tubes

Panel 2: 20-plex, FWD and REV pools provided as preformulated oligos in tubes

Final pooled panel size: 270 primers in each forward and reverse pool

- a. Determine the individual primer concentrations needed for each panel (See [Table 1](#) or pools with <500 primers):
 - i. Panel 1: 100 nM per primer.
 - ii. Panel 2: 250 nM per primer.
- b. Dilute Panel 2 to 100 nM per primer:
Dilute stock by combining 50 μL of Panel 2 with 75 μL IDTE, pH 7.5.
- c. Combine equal volumes for each FWD and REV pool (final panel = 50 nM/primer):
 - i. Combine 50 μL Panel 1 FWD + 50 μL diluted Panel 2 FWD.
 - ii. Combine 50 μL Panel 1 REV + 50 μL diluted Panel 2 REV.
- d. Proceed directly to "Perform Targeted rhAmpSeq PCR 1" in your protocol [1] or [2].

Example 2

Panel 1: 1000-plex, FWD and REV pools provided at 10X concentration

Panel 2: 2500-plex, FWD and REV pools provided at 10X concentration

Final pooled panel size: 3500 primers in each forward and reverse pools

- a. Calculate individual primer concentration for each panel:
 - i. Panel 1: 50 $\mu\text{M}/1000$ primers = 50 nM per primer.
 - ii. Panel 2: 50 $\mu\text{M}/2500$ primers = 20 nM per primer.
- b. Dilute Panel 1 to 20 nM per primer:
Dilute stock by combining 50 μL of Pool 1 with 75 μL of IDTE, pH 7.5.
- c. Combine equal volumes for each FWD and REV pool (final panel = 10 nM/primer):
 - i. Combine 50 μL diluted Panel 1 FWD + 50 μL Panel 2 FWD.
 - ii. Combine 50 μL diluted Pool 1 REV + 50 μL Panel 2 REV.
- d. Proceed directly to "Perform Targeted rhAmpSeq PCR 1" in your protocol [1] or [2].

Example 3

Panel 1: 1000-plex, FWD and REV pools provided at 10X concentration

Panel 2: 200-plex, FWD and REV pools provided at 10X concentration

Final pooled panel size: 1200 primers in each forward and reverse pool

- a. Determine the individual primer concentrations needed for each panel (See [Table 1](#) or pools with <500 primers):
 - i. Panel 1: 50 $\mu\text{M}/1000$ primers = 50 nM per primer.
 - ii. Panel 2: 100 nM per primer.
- b. Dilute Panel 2 to 50 nM per primer:
Dilute stock by combining 50 μL of Panel 2 with 50 μL of IDTE, pH 7.5.
- c. Combine equal volumes for each FWD and REV pool (final pooled panel = 25 nM/primer):
 - i. Combine 50 μL Panel 1 FWD + 50 μL diluted Panel 2 FWD.
 - ii. Combine 50 μL Panel 1 REV + 50 μL diluted Panel 2 REV.
- d. Proceed directly to "Perform Targeted rhAmpSeq PCR 1" in your protocol [1] or [2].

Scenario 3: Combining preformulated oligos in tubes and individual oligos in plates

1. For the panel shipped in plates, resuspend rhAmpSeq oligos to 50 μM (see Scenario 1, steps 1 and 2), then make separate forward and reverse pools by combining equal volumes of each primer from all panels.
2. Calculate the individual per primer concentration in each panel.

For preformulated oligos in tubes

- a. ≥ 500 -plex panels: divide 50 μM by the number of primers in the pool (i.e., 50 $\mu\text{M}/1000$ primers = 50 nM/primer).
- b. ≤ 499 -plex panels: use [Table 1](#) to determine 10X per primer concentration of each panel.

For 50 μM oligos in plates

- a. Divide 50 μM by the number of primers in the pool (i.e., 50 $\mu\text{M}/1000$ primers = 50 $\mu\text{M}/\text{primer}$).
3. Determine which of the panels has the lower per primer concentration, and dilute other panels to that concentration.
4. Combine equal volumes of each panel to create a single pooled panel, then proceed directly to "Perform Targeted rhAmpSeq PCR 1" in your protocol [1] or [2].



Note: Resulting per primer concentration will be lower than amount recommended in [Table 1](#); reaction and thermal cycling parameters may require optimization for optimal assay performance.

Example 1

Panel 1: 80-plex, 0.4 nmol scale shipped dried down in plates

Panel 2: 1000-plex, forward (FWD) and reverse (REV) pools provided lab ready at 10X concentration in tubes

Final pooled panel size: 1080 primers in each forward and reverse pool

- a. For Panel 1, resuspend each primer to 50 μM and create separate forward and reverse pools:
 - i. Add 8 μL of IDTE, pH 7.5 buffer to each well.
 - ii. Combine 2 μL of each forward primer in the FWD pool and 2 μL of each reverse primer in the REV pool (total of 160 μL per pool).

- b. Calculate individual primer concentration for each panel.
 - i. Panel 1: 50 $\mu\text{M}/80$ primers = 625 nM per primer.
 - ii. Panel 2: 50 $\mu\text{M}/1000$ primers = 50 nM per primer.
- c. Dilute Panel 1 to 50 nM per primer:
Dilute stock by combining 50 μL of Panel 1 with 575 μL of IDTE, pH 7.5.
- d. Combine equal volumes for each FWD and REV pool (final panel = 25 nM/primer):
 - i. Combine 50 μL diluted Panel 1 FWD + 50 μL Panel 2 FWD.
 - ii. Combine 50 μL diluted Panel 1 REV + 50 μL Panel 2 REV.
- e. Proceed directly to "Perform Targeted rhAmpSeq PCR 1" in your protocol [1] or [2].

References

1. Integrated DNA Technologies. (2019) **Protocol: rhAmpSeq Library Preparation for Targeted Amplicon Sequencing.** [Online] Coralville, IA, Integrated DNA Technologies, Inc. [Accessed 21 April, 2020]
2. Integrated DNA Technologies. (2019) **Protocol: High-throughput rhAmpSeq Library Preparation for Targeted Amplicon Sequencing.** [Online] Coralville, IA, Integrated DNA Technologies, Inc. [Accessed 21 April, 2020]

Technical support: applicationsupport@idtdna.com

For more than 30 years, IDT's innovative tools and solutions for genomics applications have been driving advances that inspire scientists to dream big and achieve their next breakthroughs. IDT develops, manufactures, and markets nucleic acid products that support the life sciences industry in the areas of academic and commercial research, agriculture, medical diagnostics, and pharmaceutical development. We have a global reach with personalized customer service.

> SEE WHAT MORE WE CAN DO FOR YOU AT WWW.IDTDNA.COM.

For Research Use Only.

© 2020 Integrated DNA Technologies, Inc. All rights reserved. rhAmpSeq is a trademark of Integrated DNA Technologies, Inc., and registered in the USA. All other marks are the property of their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks. NGS-10196-QR 4/2020