


Alt-R™ HDR ENHANCER V2 IMPROVES EFFICIENCY OF HOMOLOGY-DIRECTED REPAIR (HDR) IN CELLS TRANSFECTED BY LIPOFECTION

CRISPR-mediated HDR efficiency varies by cell line, editing site, and the desired insert. IDT's internal experimental results demonstrated the versatility of the Alt-R HDR Enhancer V2* to increase HDR efficiency in multiple cell lines when using electroporation to deliver different Cas enzymes with their guide RNAs. The following data shows how Alt-R HDR Enhancer V2 also showed improved HDR efficiency when cells are transfected by lipofection.



Alt-R HDR Enhancer V2 improves HDR efficiency in cells transfected by lipofection. HEK-293 cells stably expressing Cas9 were reverse transfected with gRNA complex (Alt-R CRISPR-Cas9 crRNA complexed with tracrRNA) targeting SAA1, STAT3, SERPINC1, and HPRT 38087 in the human genome (final concentration = 10 nM) and Alt-R HDR Donor Blocks designed to insert six bases at the Cas9 cleavage site (final concentration = 3 nM) using 0.75 µL of Lipofectamine® RNAiMAX® reagent (Thermo Fisher Scientific). Immediately after lipofection, cells were cultured in media containing either no treatment (dark blue), DMSO (vehicle control, light blue), 30 µM Alt-R HDR Enhancer (dark gray), 0.5 µM Alt-R HDR Enhancer V2 (light gray), or 1 µM Alt-R HDR Enhancer V2 (green). After 24 hours, the old culture media was removed and replaced with fresh cell culture media without DMSO, Alt-R HDR Enhancer, or Alt-R HDR Enhancer V2. Genomic DNA was isolated from each sample type 48 hours after lipofection, and the targeted editing sites were amplified by PCR. The efficiency was quantified using the rhAmpSeq™ CRISPR Analysis System that uses NGS analysis to assess the percentage of HDR for on-target locations. NGS reads were analyzed using the rhAmpSeq CRISPR Analysis Tool. The highest HDR was achieved using the Alt-R HDR Enhancer V2 with a final concentration of 1 µM.

* For research use only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

©2021 Integrated DNA Technologies, Inc. All rights reserved. Alt-R and rhAmpSeq are trademarks of Integrated DNA Technologies, Inc. and are registered in the US. MiSeq is a registered trademark of Illumina, Inc. Lipofectamine and RNAiMAX are registered trademarks of Thermo Fisher Scientific. All other marks are the property of their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks. MAPPs Doc ID# RUO21-0216_001 08/21

> SEE WHAT MORE WE CAN DO FOR YOU AT WWW.IDTDNA.COM.



custom oligos • qPCR • next generation sequencing • RNAi • genes & gene fragments • CRISPR genome editing