

FusionPlex-HT Lymphoma

Description

The FusionPlex-HT Lymphoma panel is a balanced pool of gene-specific primer (GSP) oligonucleotides that is optimized for use with FusionPlex-HT reagents and molecular barcode (MBC) adapters to produce targeted NGS libraries. This product insert should be used in conjunction with FusionPlex-HT protocol for Illumina® (RA-DOC-049).

FusionPlex-HT Lymphoma contains **716** GSPs targeting **125** genes commonly mutated in hematological malignancies.

Description	Part number	Storage
FusionPlex-HT Lymphoma GSP1, 24 reactions or FusionPlex-HT Lymphoma GSP1, 96 reactions	SA0077241 or SA0077961	-20°C ± 10°C
FusionPlex-HT Lymphoma GSP2, 24 reactions or FusionPlex-HT Lymphoma GSP2, 96 reactions	SA0077242 or SA0077962	
10X VCP Primer Mix, 24 reactions or 10X VCP Primer Mix, 96 reactions	SA0840 or SA0841	

Recommended PCR cycling

	Step	Temperature (°C)	Time	Cycles
First PCR reaction	1	95	3 min	1
	2	95	30 sec	15
	3	65	5 min (100% ramp rate)	
	4	72	3 min	
	5	4	Hold	1

Recommended PCR cycling (cont.)

	Step	Temperature (°C)	Time	Cycles
Second PCR reaction	1	95	3 min	1
	2	95	30 sec	20†
	3	65	5 min (100% ramp rate)	
	4	72	3 min	1
	5	4	Hold	1

†The number of PCR2 cycles may be decreased if you regularly experience library yields greater than 200 nM.

Recommended reads and multiplexing

FusionPlex-HT Lymphoma libraries should be sequenced to a minimum of **2M** reads. Starting read depth recommendations for standard profiling may be adjusted based on user needs.

Archer™ Analysis settings

Sequencing data should be processed using Archer Analysis (v7, or greater). The FusionPlex-HT Lymphoma panel requires selection of the **Fusion** pipeline and (optional) **SNV/Indel** pipeline, found under the **RNA** Input Type. See the Archer Analysis User Guide for more details on setting up your analysis.

Processing of FusionPlex-HT Lymphoma libraries requires a one-time upload of the Custom Panel GTF. Files can be obtained by contacting archer-tech@idtdna.com.

Assay targets

Gene	Accession	Exon	Variant Type	Description*
<i>AICDA</i>	NM_020661	N/A	Expression	N/A
<i>AKT3</i>	NM_005465	N/A	Mutation	K172
<i>ALK</i>	NM_004304	N/A	Mutation	T1151-C1156, F1174, L1196-S1206, G1269
<i>ALK</i>	NM_004304	2, 4, 6, 10, 16, 17, 18, 19, 20, 21, 22, 23	Fusion	5'

Gene	Accession	Exon	Variant Type	Description*
<i>ASB13</i>	NM_024701	N/A	Expression	N/A
<i>BATF3</i>	NM_018664	N/A	Expression	N/A
<i>BAX</i>	NM_138761	N/A	Mutation	G179
<i>BCL2</i>	NM_000633	N/A	Mutation	F104
<i>BCL2</i>	NM_000633	N/A	Expression	N/A
<i>BCL2</i>	NM_000633	3	Fusion	3'
<i>BCL2</i>	NM_000633	2	Fusion	5'
<i>BCL2A1</i>	NM_004049	N/A	Expression	N/A
<i>BCL3</i>	NM_005178	N/A	Expression	N/A
<i>BCL6</i>	NM_001706	N/A	Expression	N/A
<i>BCL6</i>	NM_001706	2, 3	Fusion	5'
<i>BCR</i>	NM_004327	1, 2, 3, 8, 12, 13, 14, 15, 16	Fusion	3'
<i>BIRC3</i>	NM_001165	N/A	Mutation	Q547
<i>BIRC3</i>	NM_001165	4, 5, 6, 7	Fusion	3'
<i>BLNK</i>	NM_013314	N/A	Expression	N/A
<i>BMF</i>	NM_033503	N/A	Expression	N/A
<i>BMP7</i>	NM_001719	N/A	Expression	N/A
<i>BRAF</i>	NM_004333	N/A	Mutation	V600
<i>BTK</i>	NM_000061	N/A	Mutation	C481

Gene	Accession	Exon	Variant Type	Description*
<i>CARD11</i>	NM_032415	N/A	Mutation	D230
<i>CBFB</i>	NM_022845	4,5	Fusion	3'
<i>CCDC50</i>	NM_174908	N/A	Expression	N/A
<i>CCND1</i>	NM_053056	N/A	Expression	N/A
<i>CCND1</i>	NM_053056	5	Fusion	3'
<i>CCND1</i>	NM_053056	N/A	Mutation	E36, V42-C47
<i>CCND2</i>	NM_001759	N/A	Expression	N/A
<i>CCND3</i>	NM_001760	N/A	Expression	N/A
<i>CCND3</i>	NM_001760	2	Fusion	5'
<i>CD274</i>	NM_014143	N/A	Expression	N/A
<i>CD44</i>	NM_000610	N/A	Expression	N/A
<i>CD79B</i>	NM_000626	N/A	Mutation	Y196
<i>CDC25A</i>	NM_001789	N/A	Expression	N/A
<i>CDK6</i>	NM_001259	1, 2, 3, 4	Fusion	3'
<i>CDKN2A</i>	NM_000077	N/A	Expression	N/A
<i>CDKN2B</i>	NM_004936	N/A	Expression	N/A
<i>CEBPD</i>	NM_005195	N/A	Expression	N/A
<i>CEBPE</i>	NM_001805	N/A	Expression	N/A
<i>CEBPG</i>	NM_001806	N/A	Expression	N/A

Gene	Accession	Exon	Variant Type	Description*
<i>CHIC2</i>	NM_012110	1, 2, 3	Fusion	3'
<i>CIITA</i>	NM_000246	1, 2	Fusion	3'
<i>CREB3L2</i>	NM_194071	N/A	Expression	N/A
<i>CREBBP</i>	NM_004380	N/A	Mutation	P1053, C1240, R1446, S1680-L1681
<i>CREBBP</i>	NM_004380	2, 3, 4, 5, 6	Fusion	5'
<i>CTLA4</i>	NM_005214	N/A	Expression	N/A
<i>CYB5R2</i>	NM_016229	N/A	Expression	N/A
<i>DEK</i>	NM_003472	2, 3	Fusion	3'
<i>DENND3</i>	NM_014957	N/A	Expression	N/A
<i>DLEU1</i>	NR_002605	N/A	Expression	N/A
<i>DNMT3A</i>	NM_175629	N/A	Mutation	R688-K693, R882
<i>DNMT3B</i>	NM_006892	N/A	Expression	N/A
<i>DNTT</i>	NM_004088	N/A	Expression	N/A
<i>DUSP22</i>	NM_020185	1, 2	Fusion	3'
<i>E2F2</i>	NM_004091	N/A	Expression	N/A
<i>EIF4A1</i>	NM_001416	2, 3	Fusion	5'
<i>ENTPD1</i>	NM_001776	N/A	Expression	N/A
<i>ETV6</i>	NM_001987	N/A	Mutation	Y104-R105
<i>ETV6</i>	NM_001987	1, 2, 3, 4, 5, 6	Fusion	3'

Gene	Accession	Exon	Variant Type	Description*
<i>ETV6</i>	NM_001987	2, 3, 4, 5, 6	Fusion	5'
<i>EXOC2</i>	NM_018303	N/A	Expression	N/A
<i>EZH2</i>	NM_004456	N/A	Mutation	Y646, R690
<i>FAM216A</i>	NM_013300	N/A	Expression	N/A
<i>FBXW7</i>	NM_033632	N/A	Mutation	R224-T226, R338-I347, T385, R465, R479, R505
<i>FGFR1</i>	NM_023110	12, 17	Fusion	3'
<i>FGFR1</i>	NM_023110	2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 17	Fusion	5'
<i>FOXP1</i>	NM_032682	N/A	Expression	N/A
<i>FUT8</i>	NM_004480	N/A	Expression	N/A
<i>IDH1</i>	NM_005896	N/A	Mutation	R132
<i>IDH2</i>	NM_002168	N/A	Mutation	R140, R172
<i>IL16</i>	NM_004513	N/A	Expression	N/A
<i>IRF4</i>	NM_002460	N/A	Expression	N/A
<i>IRF8</i>	NM_002163	N/A	Expression	N/A
<i>ITPKB</i>	NM_002221	N/A	Expression	N/A
<i>JAK1</i>	NM_002227	N/A	Mutation	V658, S703, R724
<i>JAK2</i>	NM_004972	N/A	Mutation	F537-F547, V617-C618, L681-R683, L855, V863, A880, V911, M929- R938, I960, R980-E985, D994
<i>JAK2</i>	NM_004972	6, 7, 8, 9, 10, 11, 12, 13, 15, 16, 17, 18, 19, 20	Fusion	5'

Gene	Accession	Exon	Variant Type	Description*
<i>JAK2</i>	NM_004972	9, 10, 11, 12	Fusion	3'
<i>JAK3</i>	NM_000215	N/A	Mutation	M511, A572-A573, R657, S789
<i>KIAA0101</i>	NM_014736	N/A	Expression	N/A
<i>KMT2A</i>	NM_005933	2, 3	Fusion	5'
<i>KMT2A</i>	NM_005933	4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35	Fusion	3'
<i>KRAS</i>	NM_004985	N/A	Mutation	G12-G13, Q61, A146
<i>LIMD1</i>	NM_014240	N/A	Expression	N/A
<i>LMO2</i>	NM_005574	N/A	Expression	N/A
<i>LRMP</i>	NM_006152	N/A	Expression	N/A
<i>LZTS1</i>	NM_021020	N/A	Expression	N/A
<i>MAL</i>	NM_002371	N/A	Expression	N/A
<i>MALT1</i>	NM_006785	9	Fusion	3'
<i>MAML3</i>	NM_018717	N/A	Expression	N/A
<i>MKL1</i>	NM_020831	4, 5, 6	Fusion	5'
<i>MLF1</i>	NM_022443	2, 3, 4	Fusion	5'
<i>MLLT10</i>	NM_004641	2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18	Fusion	5'
<i>MLLT10</i>	NM_004641	7, 8, 9, 10	Fusion	3'
<i>MME</i>	NM_000902	N/A	Expression	N/A

Gene	Accession	Exon	Variant Type	Description*
<i>MUC1</i>	NM_002456	N/A	Expression	N/A
<i>MYBL1</i>	NM_001080416	N/A	Expression	N/A
<i>MYC</i>	NM_002467	N/A	Expression	N/A
<i>MYC</i>	NM_002467	1, 2	Fusion	5'
<i>MYD88</i>	NM_002468	N/A	Mutation	V217-S219, M232, S243, L265
<i>NEK6</i>	NM_014397	N/A	Expression	N/A
<i>NFKB1</i>	NM_003998	N/A	Expression	N/A
<i>NFKB2</i>	NM_002502	14, 15, 16, 17, 18, 19, 20, 21	Fusion	3'
<i>NME1</i>	NM_000269	N/A	Expression	N/A
<i>NOTCH1</i>	NM_017617	N/A	Mutation	L1574, V1578, L1585, F1592-L1593, R1598-L1600, L1678-I1680, P2514-E2515, P2525
<i>NOTCH1</i>	NM_017617	24, 25, 26, 27, 28, 29	Fusion	5'
<i>NOTCH1</i>	NM_017617	24	Fusion	3'
<i>NOTCH1</i>	NM_017617	34	Exon Skipping	N/A
<i>NOTCH2</i>	NM_024408	N/A	Mutation	R2400
<i>NRAS</i>	NM_002524	N/A	Mutation	G12-G13, G60-Q61
<i>P2RY8</i>	NM_178129	1	Fusion	3'
<i>PAICS</i>	NM_006452	N/A	Expression	N/A
<i>PDCD1</i>	NM_005018	N/A	Expression	N/A
<i>PDCD1LG2</i>	NM_025239	1, 2, 3	Fusion	5'

Gene	Accession	Exon	Variant Type	Description*
<i>PDCD1LG2</i>	NM_025239	5, 6	Fusion	3'
<i>PDCD1LG2</i>	NM_025239	N/A	Expression	N/A
<i>PDGFRA</i>	NM_006206	N/A	Mutation	T674
<i>PDGFRA</i>	NM_006206	9, 10, 11, 12, 13, 14	Fusion	5'
<i>PIM1</i>	NM_002648	N/A	Expression	N/A
<i>PIM2</i>	NM_006875	N/A	Expression	N/A
<i>PLCG1</i>	NM_002660	N/A	Mutation	S345
<i>PLCG2</i>	NM_002661	N/A	Mutation	R665, L845
<i>PPAT</i>	NM_002703	N/A	Expression	N/A
<i>PRDM16</i>	NM_022114	1, 2, 3, 4	Fusion	5'
<i>PRKAR2B</i>	NM_002736	N/A	Expression	N/A
<i>PTPN1</i>	NM_002827	N/A	Expression	N/A
<i>PYCR1</i>	NM_006907	N/A	Expression	N/A
<i>RAB29</i>	NM_003929	N/A	Expression	N/A
<i>RAG1</i>	NM_000448	N/A	Expression	N/A
<i>RAG2</i>	NM_000536	N/A	Expression	N/A
<i>RANBP1</i>	NM_002882	N/A	Expression	N/A
<i>RHOA</i>	NM_001664	N/A	Mutation	G17
<i>S1PR2</i>	NM_004230	N/A	Expression	N/A

Gene	Accession	Exon	Variant Type	Description*
<i>SERPINA9</i>	NM_175739	N/A	Expression	N/A
<i>SF3B1</i>	NM_012433	N/A	Mutation	E622-N626, H662-K666, K700, G742
<i>SH3BP5</i>	NM_004844	N/A	Expression	N/A
<i>STAT3</i>	NM_003150	N/A	Mutation	Y640-N647, D661
<i>STAT5B</i>	NM_012448	N/A	Mutation	N642
<i>STAT6</i>	NM_003153	N/A	Mutation	D419
<i>STIL</i>	NM_003035	1, 2	Fusion	3'
<i>STRBP</i>	NM_018387	N/A	Expression	N/A
<i>TCF3</i>	NM_003200	11, 12, 13, 14, 15, 16, 17, 18	Fusion	3'
<i>TNFRSF13B</i>	NM_012452	N/A	Expression	N/A
<i>TNFSF4</i>	NM_003326	N/A	Expression	N/A
<i>TP63</i>	NM_003722	3, 4, 5	Fusion	5'
<i>WT1</i>	NM_000378	N/A	Mutation	Q301-T303, G352-Q354, V359-T365, Q445-F448
<i>WT1</i>	NM_000378	N/A	Expression	N/A
<i>XPO1</i>	NM_003400	N/A	Mutation	E571

*The mutations listed under the Description column are targeted by the assay design. Version 6.2 and earlier of Archer Analysis may not support RNA SNV/Indel variant calling at exon junctions depending on the sequence context (SNVs ≤5bp, indels ≤30bp). *De Novo* RNA SNV/Indel and Internal Tandem Duplication mutation detection are not supported on the Ion Torrent Sequencing System.

SNPs and sites targeted for sample tracking

rs560681	rs430046	rs987640	rs10776839	rs12393891
rs740598	rs8078417	rs6444724	rs6530357	chrX 4429309
rs1498553	rs9951171	rs6811238	rs5971553	chrX 11314433



Product Insert

FusionPlex™-HT Lymphoma panel

rs10773760	rs576261	rs13182883	rs5953060	chrY 6738552
rs1058083	rs1109037	rs214955	rs6524626	chrY 19490214
rs4530059	rs1523537	rs321198	rs5940270	
rs1821380	rs221956	rs4606077	rs722847	

SNPs may be used in combination to uniquely tag and track samples over time. Contact archer-tech@idtdna.com for further details.

Limitations of use

For research use only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

Safety data sheets pertaining to this product are available upon request.

This product or service is licensed under one or more of the following U.S. Patents: 8,835,358; 9,290,808; 9,290,809; 9,315,857; 9,708,659; and 9,816,137 owned by BD and is licensed solely for the use described in the associated product literature. No other rights, implied or otherwise, are granted to purchaser hereunder. Purchaser agrees, by way of example and not limitation, not to use this product to trace back the origin of a nucleic acid to an individual cell as a discrete entity (e.g., single cell analysis).

© 2023 Integrated DNA Technologies, Inc. All rights reserved. FusionPlex, VariantPlex, LiquidPlex, Immunoverse, Archer Analysis, and Archer Assay Marketplace are trademarks of Integrated DNA Technologies, Inc. All other marks are the property of their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks.