

VariantPlex-HT Core Myeloid

Description

The VariantPlex-HT Core Myeloid panel is a balanced pool of gene-specific primer (GSP) oligonucleotides that is optimized for use with VariantPlex-HT reagents and molecular barcode (MBC) adapters to produce targeted NGS libraries. This product insert should be used in conjunction with VariantPlex-HT protocol for Illumina® (RA-DOC-058).

VariantPlex-HT Core Myeloid contains **748** GSPs targeting **37** genes commonly mutated in myeloid malignancies.

Description	Part number	Storage
VariantPlex-HT Core Myeloid GSP1, 24 reactions or VariantPlex-HT Core Myeloid GSP1, 96 reactions	SA5030241 or SA5030961	-20°C ± 10°C
VariantPlex-HT Core Myeloid GSP2, 24 reactions or VariantPlex-HT Core Myeloid GSP2, 96 reactions	SA5030242 or SA5030962	

Required reagent volumes

Protocol reference	Protocol step	Reagent	Volume per reaction (µL)
A	Ligation Step 2 Elution	5mM NaOH	20
B	First PCR	VariantPlex-HT Core Myeloid GSP1	8
C	First PCR	10mM Tris-HCl pH 8.0	18
D	First PCR	Purified PCR1 eluate	16
E	Second PCR	VariantPlex-HT Core Myeloid GSP2	8

Recommended PCR cycling

	Step	Temperature (°C)	Time	Cycles
First PCR reaction	1	95	3 min	1
	2	95	30 sec	16
	3	63	5 min (100% ramp rate)	
	4	72	3 min	1
	5	4	Hold	1
Second PCR reaction	1	95	3 min	1
	2	95	30 sec	20 [†]
	3	65	5 min (100% ramp rate)	
	4	72	3 min	1
	5	4	Hold	1

†The number of PCR2 cycles may be decreased if you regularly experience library yields greater than 200 nM.

Recommended reads and multiplexing

VariantPlex-HT Core Myeloid libraries should be sequenced to a minimum of **3M** reads. Starting read depth recommendations for standard profiling may be adjusted based on user needs.

Archer™ Analysis settings

Sequencing data should be processed using Archer Analysis (v7, or greater). The VariantPlex-HT Core Myeloid panel requires selection of the **SNV/Indel, Copy Number Variation, and Structural Variation** pipelines, found under the **DNA** Input Type (see the Archer Analysis User Guide for more details on setting up your analysis). Selection of the DNA Target Coverage pipeline is optional.

Processing of VariantPlex-HT Core Myeloid libraries requires a one-time upload of the Panel GTF. When performing DNA Target Coverage analysis, users must also select a Region of Interest BED file. Users may optionally add a Targeted Mutations VCF file for targeted SNV/Indel detection. Files can be obtained by contacting archer-tech@idtdna.com

Assay targets

Gene	Accession	Exon
<i>ABL1</i>	NM_005157	4,5,6,7,8,9,10
<i>ANKRD26</i>	NM_014915	1 (c.-113-c.-134)
<i>ASXL1</i>	NM_015338	11,12,13
<i>BCOR</i>	NM_017745	2,3,4,5,6,7,9,10,11,12,13,14,15
<i>BCOR</i>	NM_001123385	8
<i>BRAF</i>	NM_004333	11,15
<i>CALR</i>	NM_004343	8,9
<i>CBL</i>	NM_005188	8,9
<i>CEBPA</i>	NM_004364	1
<i>CSF3R</i>	NM_000760	10,14,15,16
<i>CSF3R</i>	NM_156039	17
<i>CSF3R</i>	NM_172313	18
<i>DDX41</i>	NM_016222	1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17
<i>DNMT3A</i>	NM_022552	2,3,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19,20,21,22,23
<i>DNMT3A</i>	NM_153759	1,2
<i>DNMT3A</i>	NM_175630	4
<i>ETNK1</i>	NM_018638	3
<i>ETV6</i>	NM_001987	1,2,3,4,5,6,7,8
<i>EZH2</i>	NM_004456	2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19,20
<i>FLT3</i>	NM_004119	8,9,10,11,12,13,14,15,16,17,19,20,21
<i>GATA1</i>	NM_002049	2
<i>GATA2</i>	NM_032638	2,3,4,5,6
<i>IDH1</i>	NM_005896	3,4

Gene	Accession	Exon
<i>IDH2</i>	NM_002168	4,6
<i>JAK2</i>	NM_004972	12,13,14,15,16
<i>KIT</i>	NM_000222	2,8,9,10,11,12,13,14,15,17,18
<i>KRAS</i>	NM_004985	2,3,4
<i>MPL</i>	NM_005373	10,12
<i>NPM1</i>	NM_002520	11
<i>NRAS</i>	NM_002524	2,3,4
<i>PHF6</i>	NM_032458	9,10
<i>PHF6</i>	NM_032335	2,3,4,5,6,7,8
<i>PTPN11</i>	NM_002834	3,4,7,8,12,13
<i>PTPN11</i>	NM_080601	11
<i>RUNX1</i>	NM_001754	2,3,5,6,7,8,9
<i>RUNX1</i>	NM_001122607	1,5
<i>SETBP1</i>	NM_015559	4 (p.799-p.950)
<i>SF3B1</i>	NM_012433	13,14,15,16,17,18
<i>SRSF2</i>	NM_003016	1,2
<i>STAG2</i>	NM_006603	2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19,20,21,22,23,24,25,26,27,28,29,30,31,32,33
<i>STAG2</i>	NM_001042749	32
<i>TET2</i>	NM_001127208	4,5,6,7,8,9,10,11
<i>TET2</i>	NM_017628	3
<i>TP53</i>	NM_000546	1,2,3,4,5,6,7,8,9,10,11
<i>TP53</i>	NM_001276696	10
<i>TP53</i>	NM_001276695	10
<i>U2AF1</i>	NM_006758	2,5,6

Gene	Accession	Exon
WT1	NM_000378	1,2,3,4,5,6,7,9
WT1	NM_001198552	8
ZRSR2	NM_005089	1,2,3,4,5,6,7,8,9,10,11

Genes targeted for CNV

ASXL1	ETV6	RUNX1	U2AF1	BCOR	EZH2
TET2	WT1	CBL	FLT3	TP53	ZRSR2

Please contact archer-tech@idtdna.com to inquire about enabling additional genes for CNV detection.

SNPs and sites targeted for sample tracking

rs560681	rs430046	rs987640	rs10776839	rs12393891
rs740598	rs8078417	rs6444724	rs6530357	chrX 4429309
rs1498553	rs9951171	rs6811238	rs5971553	chrX 11314433
rs10773760	rs576261	rs13182883	rs5953060	chrY 6738552
rs1058083	rs1109037	rs214955	rs6524626	chrY 19490214
rs4530059	rs1523537	rs321198	rs5940270	
rs1821380	rs221956	rs4606077	rs722847	

SNPs may be used in combination to uniquely tag and track samples over time. Contact archer-tech@idtdna.com for further details.



Product Insert

VariantPlex™-HT Core Myeloid panel

Limitations of use

For research use only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

Safety data sheets pertaining to this product are available upon request.

This product or service is licensed under one or more of the following U.S. Patents: 8,835,358; 9,290,808; 9,290,809; 9,315,857; 9,708,659; and 9,816,137 owned by BD and is licensed solely for the use described in the associated product literature. No other rights, implied or otherwise, are granted to purchaser hereunder. Purchaser agrees, by way of example and not limitation, not to use this product to trace back the origin of a nucleic acid to an individual cell as a discrete entity (e.g., single cell analysis).

© 2023 Integrated DNA Technologies, Inc. All rights reserved. FusionPlex, VariantPlex, LiquidPlex, Immunoverse, Archer Analysis, and Archer Assay Marketplace are trademarks of Integrated DNA Technologies, Inc. All other marks are the property of their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks.