

28 April 2020

Dear Valued Customer:

This letter serves as notification of a change to the 2019-nCoV CDC qPCR Probe Assay kit manufactured by IDT for use under the CDC Emergency Use Authorization (EUA).

As part of IDT's quality management system responsibility to its valued customers you are hereby notified of a change to the 2019-nCoV CDC qPCR Probe Assay kits manufactured by IDT for use under the CDC EUA. Pursuant to an amendment to the CDC Emergency Use Authorization dated March 8th, 2020; the Iowa Black-ZEN technology was approved for use in the CDC protocol as an equivalent alternative to the quenchers used in the original CDC test panel publication.

Going forward, CDC EUA kits provided for sale by IDT may include the use of the Iowa Black-ZEN technology. EUA CDC lots produced with the Iowa Black-ZEN quencher probe technology will still follow all protocols for testing and release by CDC prior to being made available for sale.

For any concerns you may have related to this change, please contact Scientific Applications Support at applicationsupport@idtdna.com.

Best Regards,



Amy Stanlick
Director Global Quality Assurance

www.idtdna.com

Integrated DNA Technologies, Inc.
1710 Commercial Park
Coralville, Iowa 52241
USA

